

Arthrogryposes Panel

Test code: MA0501

Is a 78 gene panel that includes assessment of non-coding variants.

Is ideal for patients with a clinical suspicion of arthrogryposis or fetal akinesia.

Arthrogryposes are a group of disorders that involve congenital joint contractures. This comprehensive panel includes Fetal Akinesia Deformation Sequence / LMPS / Related Disorder Panel and covers, but is not limited to the disorders covered by the subpanels. This panel enables effective differential diagnostics of arthrogryposes and associated diseases. The Panel is also part of Comprehensive Skeletal / Malformation Syndrome Panel.

About Arthrogryposes

Arthrogryposis (also known as arthrogryposis multiplex congenita, AMC) is characterized by congenital contractures of 2 or more different body areas without a primary neurologic or muscle disease. Children born with joint contractures have abnormal fibrosis of the muscle tissue causing muscle shortening, and therefore are unable to perform passive extension and flexion in the affected joints. Arthrogryposis has been divided into three groups: amyoplasia, distal arthrogryposis, and syndromic. Amyoplasia is characterized by severe joint contractures and muscle weakness while distal arthrogryposis mainly involves the hands and feet. Syndromic arthrogryposis consists with a primary neurological or muscle disease. 70-80% of arthrogryposes are caused by neurological abnormalities and most types that have primary neurological or muscle disease result from an underlying genetic syndrome. More than 35 specific genetic disorders associated with arthrogryposis have been described. Fetal akinesia deformation sequence syndrome (FADS) is characterised by decreased fetal movement (fetal akinesia) as well as intrauterine growth restriction, arthrogryposis, and developmental anomalies.

Availability

Results in 3-4 weeks

Gene set description

Genes in the Arthrogryposes Panel and their clinical significance

| Gene | Associated phenotypes | Inheritance | ClinVar | HGMD |
|---------|---|-------------|---------|------|
| ACTA1 | Myopathy | AD/AR | 68 | 212 |
| ADGRG6 | Lethal congenital contracture syndrome 9 | AR | 4 | 4 |
| AGRN | Myasthenic syndrome, congenital | AR | 14 | 16 |
| BIN1 | Myopathy, centronuclear | AD/AR | 9 | 15 |
| CACNA1E | Epileptic encephalopathy | AD | 8 | 6 |
| CASK | Mental retardation and microcephaly with pontine and cerebellar hypoplasia, FG syndrome, Mental retardation | XL | 87 | 112 |
| CFL2 | Nemaline myopathy | AR | 2 | 9 |
| CHAT | Myasthenic syndrome, congenital | AR | 24 | 73 |
| CHRNA1 | Myasthenic syndrome, congenital | AD/AR | 28 | 35 |
| CHRNB1 | Myasthenic syndrome | AD/AR | 11 | 11 |

Blueprint Genetics

| | | | | |
|------------------------|--|-------|-----|-----|
| CHRND | Myasthenic syndrome | AD/AR | 18 | 26 |
| CHRNE | Myasthenic syndrome | AD/AR | 48 | 134 |
| CHRNG | Multiple pterygium syndrome, Escobar syndrome | AR | 17 | 34 |
| CHST14 | Ehlers-Danlos syndrome, musculocontractural | AR | 15 | 21 |
| CHUK | Cocoon syndrome | AR | 2 | 5 |
| CNTNAP1 | Lethal congenital contracture syndrome 7 | AR | 10 | 21 |
| COL6A2 | Epilepsy, progressive myoclonic, Bethlem myopathy, Myosclerosis, congenital, Ullrich congenital muscular dystrophy | AD/AR | 101 | 182 |
| COLQ | Myasthenic syndrome, congenital | AR | 23 | 67 |
| DHCR24 | Desmosterolosis | AR | 6 | 9 |
| DOK7 | Myasthenic syndrome, congenital | AR | 28 | 75 |
| DPAGT1 | Congenital disorder of glycosylation, Myasthenic syndrome, congenital | AR | 16 | 32 |
| ECEL1 | Arthrogryposis | AR | 25 | 31 |
| EGR2 | Neuropathy, Dejerine-Sottas disease, Charcot-Marie-Tooth disease | AD/AR | 13 | 21 |
| ERBB3 | Lethal congenital contractural syndrome 2 | AR | 11 | 4 |
| ERCC5 | Xeroderma pigmentosum, Xeroderma pigmentosum/Cockayne syndrome | AR | 21 | 54 |
| ERCC6* | Xeroderma Pigmentosum-Cockayne Syndrome, De Sanctis-Cacchione syndrome | AD/AR | 87 | 135 |
| EXOSC3 | Pontocerebellar hypoplasia | AR | 11 | 19 |
| FBN2 | Congenital contractural arachnodactyly (Beals syndrome) | AD | 50 | 97 |
| FHL1* | Myopathy with postural muscle atrophy, Emery-Dreifuss muscular dystrophy, Reducing bod myopathy | XL | 26 | 62 |
| FKBP10 | Bruck syndrome 1, Osteogenesis imperfecta, type XI | AR | 20 | 44 |
| FKTN | Muscular dystrophy-dystroglycanopathy, Dilated cardiomyopathy (DCM), Muscular dystrophy-dystroglycanopathy (limb-girdle) | AD/AR | 45 | 58 |
| FLVCR2 | Proliferative vasculopathy and hydraencephaly-hydrocephaly syndrome | AR | 9 | 17 |
| GBA* | Gaucher disease | AR | 84 | 488 |
| GBE1 | Glycogen storage disease | AR | 36 | 70 |
| GFPT1 | Myasthenic syndrome, congenital | AR | 13 | 42 |
| GLDN | Lethal congenital contracture syndrome 11 | AR | 11 | 11 |
| GLE1 | Lethal congenital contracture syndrome, Arthrogryposis, lethal, with anterior horn cell disease | AR | 7 | 17 |
| KAT6B | Ohdo syndrome, SBBYS variant, Genitopatellar syndrome | AD | 47 | 73 |

Blueprint Genetics

| | | | | |
|-------------------------|---|-------|-----|-----|
| KIAA1109 | Craniofacial dysmorphism, skeletal anomalies, and mental retardation syndrome | AR | 7 | 16 |
| KLHL40 | Nemaline myopathy | AR | 11 | 26 |
| LGI4 | Arthrogryposis multiplex congenita, neurogenic, with myelin defect | AR | 9 | 7 |
| LMNA | Heart-hand syndrome, Slovenian, Limb-girdle muscular dystrophy, Muscular dystrophy, congenital, LMNA-related, Lipodystrophy (Dunnigan), Emery-Dreifuss muscular dystrophy, Malouf syndrome, Dilated cardiomyopathy (DCM), Mandibuloacral dysplasia type A, Progeria Hutchinson-Gilford type | AD/AR | 250 | 564 |
| MPZ | Neuropathy, Roussy-Levy syndrome, Dejerine-Sottas disease, Charcot-Marie-Tooth disease | AD | 108 | 241 |
| MTM1 | Myopathy, centronuclear | XL | 158 | 301 |
| MUSK | Myasthenic syndrome, congenital | AR | 17 | 22 |
| MYBPC1 | Arthrogryposis, Lethal congenital contractural syndrome | AD/AR | 7 | 7 |
| MYH2 | Inclusion body myopathy | AD | 24 | 24 |
| MYH3 | Arthrogryposis | AD | 21 | 45 |
| MYH8 | Carney complex variant, Arthrogryposis, distal, type 7, Trismus-pseudocamptodactyly syndrome | AD | 1 | 2 |
| NALCN | Neuroaxonal neurodegeneration, infantile, with facial dysmorphism, Congenital contractures of the limbs and face, hypotonia, and developmental delay | AD/AR | 47 | 50 |
| NEB* | Nemaline myopathy | AR | 305 | 309 |
| PIZO2* | Marden-Walker syndrome, Distal arthrogryposis | AD | 30 | 28 |
| PLOD2 | Bruck syndrome, Osteogenesis imperfecta type 3 | AR | 8 | 23 |
| PMM2 | Congenital disorder of glycosylation | AR | 76 | 128 |
| PPP3CA | Epileptic encephalopathy | AD | 8 | 11 |
| RAPSN | Myasthenic syndrome, congenital | AR | 26 | 58 |
| RARS2# | Pontocerebellar hypoplasia | AR | 23 | 37 |
| RIPK4 | Popliteal pterygium syndrome, lethal type, Bartsocas-Papas syndrome | AR | 4 | 15 |
| SCO2 | Leigh syndrome, Hypertrophic cardiomyopathy (HCM), Cardioencephalomyopathy, fatal infantile, due to cytochrome c oxidase deficiency, Myopia | AR | 42 | 37 |
| SELENON | Muscular dystrophy, rigid spine, Myopathy, congenital, with fiber-disproportion | AR | 38 | 63 |
| SMN1*,# | Spinal muscular atrophy | AR | 29 | 111 |
| SMN2*,# | Spinal muscular atrophy | AD | 1 | 9 |

| | | | | |
|--------------|---|-------|----|----|
| TGFB3 | Loeys-Dietz syndrome (Reinhoff syndrome), Arrhythmogenic right ventricular dysplasia | AD | 19 | 26 |
| TK2 | Mitochondrial DNA depletion syndrome | AR | 38 | 52 |
| TNNI2 | Arthrogryposis multiplex congenita | AD | 5 | 11 |
| TNNT1 | Nemaline myopathy | AR | 6 | 8 |
| TNNT3 | Arthrogryposis, distal, type 2B | AD | 3 | 4 |
| TPM2 | CAP myopathy, Nemaline myopathy, Arthrogryposis, distal | AD | 18 | 38 |
| <u>TPM3*</u> | CAP myopathy, Nemaline myopathy, Myopathy, congenital, with fiber-disproportion | AD | 21 | 27 |
| TRPV4 | Metatropic dysplasia, Spondyloepiphyseal dysplasia Maroteaux type, Parastremmatic dwarfism, Hereditary motor and sensory neuropathy, Spondylometaphyseal dysplasia Kozlowski type, Spinal muscular atrophy, Charcot-Marie-Tooth disease, Brachyolmia (autosomal dominant type), Familial Digital arthropathy with brachydactyly | AD | 61 | 78 |
| TSEN2 | Pontocerebellar hypoplasia | AR | 8 | 5 |
| TSEN54 | Pontocerebellar hypoplasia | AR | 23 | 21 |
| UBA1 | Spinal muscular atrophy, infantile | XL | 3 | 5 |
| VIPAS39 | Arthrogryposis, renal dysfunction, and cholestasis 2 | AR | 8 | 13 |
| VPS33B | Arthrogryposis - renal dysfunction - cholestasis | AD/AR | 17 | 58 |
| VRK1 | Pontocerebellar hypoplasia | AR | 9 | 9 |
| ZBTB42 | Lethal congenital contracture syndrome | AR | 2 | 1 |
| ZC4H2 | Wieacker-Wolff syndrome | XL | 20 | 16 |

*Some regions of the gene are duplicated in the genome. [Read more.](#)

The gene has suboptimal coverage (means <90% of the gene's target nucleotides are covered at >20x with mapping quality score (MQ>20) reads), and/or the gene has exons listed under Test limitations section that are not included in the panel as they are not sufficiently covered with high quality sequence reads.

The sensitivity to detect variants may be limited in genes marked with an asterisk (*) or number sign (#)

Gene refers to the HGNC approved gene symbol; Inheritance refers to inheritance patterns such as autosomal dominant (AD), autosomal recessive (AR), mitochondrial (mi), X-linked (XL), X-linked dominant (XLD) and X-linked recessive (XLR); ClinVar refers to the number of variants in the gene classified as pathogenic or likely pathogenic in this database ([ClinVar](#)); HGMD refers to the number of variants with possible disease association in the gene listed in Human Gene Mutation Database ([HGMD](#)). The list of associated, gene specific phenotypes are generated from [CGD](#) or Mitomap databases.

Non-coding disease causing variants covered by the panel

| Gene | Genomic location HG19 | HGVS | RefSeq | RS-number |
|-------|-----------------------|-------------|-------------|-----------|
| CHRNE | Chr17:4804936 | c.501-16G>A | NM_000080.3 | |
| CHRNE | Chr17:4806452 | c.-94G>A | NM_000080.3 | |

Blueprint Genetics

| | | | | |
|---------|-----------------|--|----------------|--------------|
| CHRNE | Chr17:4806453 | c.-95G>A | NM_000080.3 | |
| CHRNE | Chr17:4806454 | c.-96C>T | NM_000080.3 | rs748144899 |
| COL6A2 | Chr21:47538492 | c.1117-35_1118dupAAAAGACGTGAGGCTGATTCTGCAAACCCCTCCAGGG | NM_001849.3 | |
| COL6A2 | Chr21:47541407 | c.1459-63G>A | NM_001849.3 | |
| ERCC5 | Chr13:103514354 | c.881-26T>G | NM_000123.3 | |
| ERCC6 | Chr10:50681659 | c.2599-26A>G | NM_000124.3 | rs4253196 |
| EXOSC3 | Chr9:37782146 | c.475-12A>G | NM_016042.3 | rs370087266 |
| FBN2 | Chr5:127670560 | c.3974-24A>C | NM_001999.3 | |
| FBN2 | Chr5:127670562 | c.3974-26T>G | NM_001999.3 | |
| FBN2 | Chr5:127671284 | c.3725-15A>G | NM_001999.3 | |
| FKTN | Chr9:108368857 | c.648-1243G>T | NM_006731.2 | |
| GBA | Chr1:155205646 | c.1225-14_1225-11delTGTCinsAGT | NM_000157.3 | |
| GBA | Chr1:155208109 | c.589-12C>G | NM_000157.3 | |
| GBA | Chr1:155211053 | c.-150A>G | NM_000157.3 | rs1232943445 |
| GBE1 | Chr3:81542964 | c.2053-3358_2053-3350delGTGTGGTGGinsTGTTTTTACATGACAGGT | NM_000158.3 | rs869320698 |
| LMNA | Chr1:156100609 | c.513+45T>G | NM_170707.3 | |
| LMNA | Chr1:156105681 | c.937-11C>G | NM_170707.3 | rs267607645 |
| LMNA | Chr1:156107037 | c.1608+14G>A | NM_170707.3 | |
| LMNA | Chr1:156107433 | c.1609-12T>G | NM_170707.3 | rs267607582 |
| MTM1 | ChrX:149767035 | c.137-19_137-16delACTT | NM_000252.2 | |
| MTM1 | ChrX:149767045 | c.137-11T>A | NM_000252.2 | |
| MTM1 | ChrX:149783032 | c.232-26_232-23delGACT | NM_000252.2 | |
| MTM1 | ChrX:149808833 | c.529-909A>G | NM_000252.2 | |
| MTM1 | ChrX:149818176 | c.868-13T>A | NM_000252.2 | |
| NEB | Chr2:152355017 | c.24220-151C>A | NM_001271208.1 | |
| NEB | Chr2:152410918 | c.19429-381_19429-379delTTTinsA | NM_001271208.1 | |
| PMM2 | Chr16:8891573 | | NM_000303.2 | |
| PMM2 | Chr16:8898599 | c.179-25A>G | NM_000303.2 | rs760689221 |
| PMM2 | Chr16:8926102 | c.640-15479C>T | NM_000303.2 | rs1258107584 |
| PMM2 | Chr16:8941558 | c.640-23A>G | NM_000303.2 | |
| RAPSN | Chr11:47469717 | c.193-15C>A | NM_005055.4 | |
| RAPSN | Chr11:47470715 | c.-199C>G | NM_005055.4 | |
| RAPSN | Chr11:47470726 | c.-210A>G | NM_005055.4 | rs786200905 |
| RARS2 | Chr6:88244587 | c.613-3927C>T | NM_020320.3 | |
| SELENON | Chr1:26143316 | c.*1107T>C | NM_020451.2 | |
| TGFB3 | Chr14:76425035 | c.*495C>T | NM_003239.2 | rs387906514 |



| | | | | |
|--------|----------------|-------------|-------------|-------------|
| TGFB3 | Chr14:76447266 | c.-30G>A | NM_003239.2 | rs770828281 |
| VPS33B | Chr15:91550814 | c.499-11G>A | NM_018668.3 | |

Test Strengths

All exons of the *GBA* gene have segmentally duplicated pseudogenes that reduce sensitivity of NGS diagnostics in general. However, Blueprint Genetics custom assay has good coverage (>20x) with high mapping rates (mapping quality >40) for 100.0% of the target regions in *GBA* gene. Our validation showed high mean coverage of 184X for the *GBA* gene. Thus, our NGS Panel is not expected to have major limitations in detecting variants in *GBA* gene although clinical validation has not been performed at large scale for Gaucher disease.

Deletion / duplication analysis (either in isolation or as part of Plus analysis including sequencing) testing can detect the copy number of *SMN1* exon 7, which is commonly used as a marker for copy number of the *SMN1* gene. In individuals identified to have homozygous *SMN1* deletions, we include reporting of *SMN2* copy number.

The strengths of this test include:

- CAP accredited laboratory
- CLIA-certified personnel performing clinical testing in a CLIA-certified laboratory
- Powerful sequencing technologies, advanced target enrichment methods and precision bioinformatics pipelines ensure superior analytical performance
- Careful construction of clinically effective and scientifically justified gene panels
- Some of the panels include the whole mitochondrial genome (please see the Panel Content section)
- Our Nucleus online portal providing transparent and easy access to quality and performance data at the patient level
- Our publicly available analytic validation demonstrating complete details of test performance
- ~2,000 non-coding disease causing variants in our clinical grade NGS assay for panels (please see 'Non-coding disease causing variants covered by this panel' in the Panel Content section)
- Our rigorous variant classification scheme
- Our systematic clinical interpretation workflow using proprietary software enabling accurate and traceable processing of NGS data
- Our comprehensive clinical statements

Test Limitations

Due to high sequence homology between *SMN1* and *SMN2*, this panel has not been validated to detect single nucleotide variants or small insertions and deletions in *SMN1* which are associated with spinal muscular atrophy in a small number of patients (<5%). The following exons are not included in the panel as they are not sufficiently covered with high quality sequence reads: *SELENON* (NM_020451:3), *TK2* (NM_001271934:3), *TSEN2* (NM_001321278:12). Genes with suboptimal coverage in our assay are marked with number sign (#) and genes with partial, or whole gene, segmental duplications in the human genome are marked with an asterisk (*) if they overlap with the UCSC pseudogene regions. Gene is considered to have suboptimal coverage when >90% of the gene's target nucleotides are not covered at >20x with mapping quality score (MQ>20) reads. The technology may have limited sensitivity to detect variants in genes marked with these symbols (please see the Panel content table above).

This test does not detect the following:

- Complex inversions
- Gene conversions
- Balanced translocations
- Some of the panels include the whole mitochondrial genome but not all (please see the Panel Content section)
- Repeat expansion disorders unless specifically mentioned
- Non-coding variants deeper than ± 20 base pairs from exon-intron boundary unless otherwise indicated (please see above Panel Content / non-coding variants covered by the panel).

This test may not reliably detect the following:

Blueprint Genetics

- Low level mosaicism in nuclear genes (variant with a minor allele fraction of 14.6% is detected with 90% probability)
- Stretches of mononucleotide repeats
- Low level heteroplasmy in mtDNA (>90% are detected at 5% level)
- Indels larger than 50bp
- Single exon deletions or duplications
- Variants within pseudogene regions/duplicated segments
- Some disease causing variants present in mtDNA are not detectable from blood, thus post-mitotic tissue such as skeletal muscle may be required for establishing molecular diagnosis.

The sensitivity of this test may be reduced if DNA is extracted by a laboratory other than Blueprint Genetics.

For additional information, please refer to the Test performance section and see our Analytic Validation.

Test performance

The Blueprint Genetics arthrogyroses panel covers classical genes associated with arthrogyrosis. The genes on the panel have been carefully selected based on scientific literature, mutation databases and our experience.

Our panels are sliced from our high-quality whole exome sequencing data. Please see our sequencing and detection performance table for different types of alterations at the whole exome level (Table).

Assays have been validated for different starting materials including EDTA-blood, isolated DNA (no FFPE), saliva and dry blood spots (filter card) and all provide high-quality results. The diagnostic yield varies substantially depending on the assay used, referring healthcare professional, hospital and country. Blueprint Genetics' Plus Analysis (Seq+Del/Dup) maximizes the chance to find a molecular genetic diagnosis for your patient although Sequence Analysis or Del/Dup Analysis may be a cost-effective first line test if your patient's phenotype is suggestive of a specific mutation type.

The genes on the panel have been carefully selected based on scientific literature, mutation databases and our experience.

Our panels are sectioned from our high-quality, clinical grade NGS assay. Please see our sequencing and detection performance table for details regarding our ability to detect different types of alterations (Table).

| | Sensitivity % (TP/(TP+FN)) | Specificity % |
|---|----------------------------|---------------|
| Single nucleotide variants | 99.89% (99,153/99,266) | >99.9999% |
| Insertions, deletions and indels by sequence analysis | | |
| 1-10 bps | 96.9% (7,563/7,806) | >99.9999% |
| 11-50 bps | 99.13% (2,524/2,546) | >99.9999% |
| Copy number variants (exon level dels/dups) | | |
| 1 exon level deletion (heterozygous) | 100% (20/20) | NA |
| 1 exon level deletion (homozygous) | 100% (5/5) | NA |
| 1 exon level deletion (het or homo) | 100% (25/25) | NA |
| 2-7 exon level deletion (het or homo) | 100% (44/44) | NA |
| 1-9 exon level duplication (het or homo) | 75% (6/8) | NA |
| Simulated CNV detection | | |
| 5 exons level deletion/duplication | 98.7% | 100.00% |

Microdeletion/-duplication sdrs (large CNVs, n=37)

Size range (0.1-47 Mb) 100% (37/37)

The performance presented above reached by Blueprint Genetics high-quality, clinical grade NGS sequencing assay with the following coverage metrics

Mean sequencing depth 143X

Nucleotides with >20x sequencing coverage (%) 99.86%

Performance of Blueprint Genetics Mitochondrial Sequencing Assay.

| | Sensitivity % | Specificity % |
|---|--------------------|---------------|
| ANALYTIC VALIDATION (NA samples; n=4) | | |
| Single nucleotide variants | | |
| Heteroplasmic (45-100%) | 100.0% (50/50) | 100.0% |
| Heteroplasmic (35-45%) | 100.0% (87/87) | 100.0% |
| Heteroplasmic (25-35%) | 100.0% (73/73) | 100.0% |
| Heteroplasmic (15-25%) | 100.0% (77/77) | 100.0% |
| Heteroplasmic (10-15%) | 100.0% (74/74) | 100.0% |
| Heteroplasmic (5-10%) | 100.0% (3/3) | 100.0% |
| Heteroplasmic (<5%) | 50.0% (2/4) | 100.0% |
| CLINICAL VALIDATION (n=76 samples) | | |
| All types | | |
| Single nucleotide variants n=2026 SNVs | | |
| Heteroplasmic (45-100%) | 100.0% (1940/1940) | 100.0% |
| Heteroplasmic (35-45%) | 100.0% (4/4) | 100.0% |
| Heteroplasmic (25-35%) | 100.0% (3/3) | 100.0% |
| Heteroplasmic (15-25%) | 100.0% (3/3) | 100.0% |
| Heteroplasmic (10-15%) | 100.0% (9/9) | 100.0% |
| Heteroplasmic (5-10%) | 92.3% (12/13) | 99.98% |
| Heteroplasmic (<5%) | 88.9% (48/54) | 99.93% |
| Insertions and deletions by sequence analysis n=40 indels | | |
| Heteroplasmic (45-100%) 1-10bp | 100.0% (32/32) | 100.0% |

Blueprint Genetics



| | | |
|---|-----------------|-------------------|
| Heteroplasmic (5-45%) 1-10bp | 100.0% (3/3) | 100.0% |
| Heteroplasmic (<5%) 1-10bp | 100.0% (5/5) | 99,997% |
| SIMULATION DATA /(mitomap mutations) | | |
| Insertions, and deletions 1-24 bps by sequence analysis; n=17 | | |
| Homoplasmic (100%) 1-24bp | 100.0% (17/17) | 99.98% |
| Heteroplasmic (50%) | 100.0% (17/17) | 99.99% |
| Heteroplasmic (25%) | 100.0% (17/17) | 100.0% |
| Heteroplasmic (20%) | 100.0% (17/17) | 100.0% |
| Heteroplasmic (15%) | 100.0% (17/17) | 100.0% |
| Heteroplasmic (10%) | 94.1% (16/17) | 100.0% |
| Heteroplasmic (5%) | 94.1% (16/17) | 100.0% |
| Copy number variants (separate artificial mutations; n=1500) | | |
| Homoplasmic (100%) 500 bp, 1kb, 5 kb | 100.0% | 100.0% |
| Heteroplasmic (50%) 500 bp, 1kb, 5 kb | 100.0% | 100.0% |
| Heteroplasmic (30%) 500 bp, 1kb, 5 kb | 100.0% | 100.0% |
| Heteroplasmic (20%) 500 bp, 1kb, 5 kb | 99.7% | 100.0% |
| Heteroplasmic (10%) 500 bp, 1kb, 5 kb | 99.0% | 100.0% |
| The performance presented above reached by following coverage metrics at assay level (n=66) | | |
| | Mean of medians | Median of medians |
| Mean sequencing depth MQ0 (clinical) | 18224X | 17366X |
| Nucleotides with >1000x MQ0 sequencing coverage (%) (clinical) | 100% | |
| rho zero cell line (=no mtDNA), mean sequencing depth | 12X | |

Bioinformatics

The target region for each gene includes coding exons and ± 20 base pairs from the exon-intron boundary. In addition, the panel includes non-coding variants if listed above (Non-coding variants covered by the panel). Some regions of the gene(s) may be removed from the panel if specifically mentioned in the "Test limitations" section above. The sequencing data generated in our laboratory is analyzed with our proprietary data analysis and annotation pipeline, integrating state-of-the-art algorithms and industry-standard software solutions. Incorporation of rigorous quality control steps throughout the workflow of the pipeline ensures the consistency, validity and accuracy of results. Our pipeline is streamlined to maximize sensitivity without sacrificing specificity. We have incorporated a number of reference population databases and mutation databases such as, but not limited, to [1000 Genomes Project](#), [gnomAD](#), [ClinVar](#) and [HGMD](#) into our clinical interpretation software to make the process effective and efficient. For missense variants, *in silico* variant prediction tools such as [SIFT](#), [PolyPhen](#), [MutationTaster](#) are used to assist with variant classification. Through our [online ordering and statement reporting system](#), [Nucleus](#), the customer has an access to details of the analysis, including patient specific sequencing metrics, a gene level coverage plot and a list of regions with inadequate coverage if present. This reflects our mission to build fully transparent



[diagnostics where customers have easy access to crucial details of the analysis process.](#)

Clinical interpretation

We provide customers with the most comprehensive clinical report available on the market. Clinical interpretation requires a fundamental understanding of clinical genetics and genetic principles. At Blueprint Genetics, our PhD molecular geneticists, medical geneticists and clinical consultants prepare the clinical statement together by evaluating the identified variants in the context of the phenotypic information provided in the requisition form. Our goal is to provide clinically meaningful statements that are understandable for all medical professionals regardless of whether they have formal training in genetics.

Variant classification is the corner stone of clinical interpretation and resulting patient management decisions. Our classifications follow the [Blueprint Genetics Variant Classification Schemes](#) based on the [ACMG guideline 2015](#). Minor modifications were made to increase reproducibility of the variant classification and improve the clinical validity of the report. Our experience with tens of thousands of clinical cases analyzed at our laboratory allowed us to further develop the industry standard.

The final step in the analysis of sequence variants is confirmation of variants classified as pathogenic or likely pathogenic using bi-directional Sanger sequencing. Variant(s) fulfilling the following criteria are not Sanger confirmed: the variant quality score is above the internal threshold for a true positive call, and visual check-up of the variant at IGV is in-line with the variant call. Reported variants of uncertain significance are confirmed with bi-directional Sanger sequencing only if the quality score is below our internally defined quality score for true positive call. Reported copy number variations with a size <10 exons are confirmed by orthogonal methods such as qPCR if the specific CNV has been seen less than three times at Blueprint Genetics.

Our clinical statement includes tables for sequencing and copy number variants that include basic variant information (genomic coordinates, HGVS nomenclature, zygosity, allele frequencies, in silico predictions, OMIM phenotypes and classification of the variant). In addition, the statement includes detailed descriptions of the variant, gene and phenotype(s) including the role of the specific gene in human disease, the mutation profile, information about the gene's variation in population cohorts and detailed information about related phenotypes. We also provide links to the references used, congress abstracts and mutation variant databases used to help our customers further evaluate the reported findings if desired. The conclusion summarizes all of the existing information and provides our rationale for the classification of the variant.

Identification of pathogenic or likely pathogenic variants in dominant disorders or their combinations in different alleles in recessive disorders are considered molecular confirmation of the clinical diagnosis. In these cases, family member testing can be used for risk stratification within the family. In the case of variants of uncertain significance (VUS), we do not recommend family member risk stratification based on the VUS result. Furthermore, in the case of VUS, we do not recommend the use of genetic information in patient management or genetic counseling.

Our interpretation team analyzes millions of variants from thousands of individuals with rare diseases. Thus, our database, and our understanding of variants and related phenotypes, is growing by leaps and bounds. Our laboratory is therefore well positioned to re-classify previously reported variants as new information becomes available. If a variant previously reported by Blueprint Genetics is re-classified, our laboratory will issue a follow-up statement to the original ordering health care provider at no additional cost.

#}

ICD codes

Commonly used ICD-10 codes when ordering the Arthrogyposes Panel

| ICD-10 | Disease |
|--------|---------------|
| Q74.3 | Arthrogyposis |

Accepted sample types

- EDTA blood, min. 1 ml
- Purified DNA, min. 3µg*
- Saliva (Oragene DNA OG-500 kit)

Label the sample tube with your patient's name, date of birth and the date of sample collection.

Note that we do not accept DNA samples isolated from formalin-fixed paraffin-embedded (FFPE) tissue.

Resources

- [American Association of Neuromuscular & Electrodiagnostic Medicine](#)
- [Arthrogryposis Multiplex Congenita Support](#)
- [Bamshad M et al. Arthrogryposis: a review and update. J Bone Joint Surg Am. 2009 Jul;91 Suppl 4:40-6.](#)
- [European Alliance for Arthrogryposis](#)
- [GeneReviews](#)
- [GeneReviews - Congenital Contractural Arachnodactyly](#)
- [GeneReviews - Congenital Myasthenic Syndromes](#)
- [GeneReviews - Spinal Muscular Atrophy](#)
- [Genetic and Rare Disease Information Center](#)
- [Kowalczyk B et al. Arthrogryposis: an update on clinical aspects, etiology, and treatment strategies. Arch Med Sci. 2016 Feb 1;12\(1\):10-24.](#)
- [Muscular Dystrophy Association](#)
- [NORD - Arthrogryposis](#)
- [NORD - Arthrogryposis Multiplex Congenita](#)
- [National Organization for Rare Disorders](#)
- [The Arthrogryposis Group](#)