

## Congenital Disorders of Glycosylation Panel

Test code: ME1901

Is ideal for patients with a clinical suspicion of a congenital disorder of N-linked glycosylation or combined defects of glycosylation affecting both the N-linked and O-linked glycosylation pathways. The genes on this panel are included in the Comprehensive Metabolism Panel.

The most types of congenital disorders of N-linked glycolysation are inherited in autosomal recessive manner. *MGAT1*, *ALG13*, *SLC35A2* and *SSR4*-related disorders are inherited in X-linked manner. In addition to congenital disorders of N-linked glycolysation, this Panel has differential diagnostics power to rare phenotypes with overlapping symptoms such as *GEN*-related myopathy and *ATP6VOA2*- related cutis laxa. This Panel is included in the Comprehensive Metabolism Panel.

### About Congenital Disorders of Glycosylation

Most subtypes of congenital disorders of glycosylation (CDG) are classified as disorders of N-glycosylation, which involves carbohydrates called N-linked oligosaccharides. These oligosaccharides are created in a specific order to create specific sugar trees, which are then attached to proteins on various cells. Disorders of N-glycosylation are due to an enzyme deficiency or other malfunction somewhere along the N-glycosylation pathway. There are 42 different enzymes in the pathway; any of them may be mutated and cause a disorder belonging to this group. Different mutated enzymes cause different phenotypes. Congenital disorders of N-linked glycosylation are a genetically and phenotypically heterogeneous group of diseases. Most commonly, symptoms begin in early infancy. Manifestations range from mild to severe, involving only protein-losing enteropathy and hypoglycemia or severe intellectual disability with malfunction of several organs. Sometimes the disorder may be fatal. Most patients require nutritional supplements. Most of the individual disorders have been observed only in a very limited number of patients. The most common ones are *PMM2*-related disorder (approximately 700 patients reported), *MPI*-related disorder (>20 patients) and *ALG6*-related disorder (>30 patients). Other types of disorder are extremely rare. In addition to congenital disorders of N-linked glycosylation, this panel has the ability to diagnose rare phenotypes with overlapping symptoms such as *GEN*-related myopathy and *ATP6VOA2*- related cutis laxa. The panel also covers genes for CDG that occur due to combined defects of glycosylation; defects affecting both the N-linked and O-linked glycosylation pathways. Genes for disorders of protein O-glycosylation, which in many cases have been classified as subtypes of other umbrella groups (e.g., muscular dystrophy), and show more dysmorphic features in general, are better known with more traditional names and can be found on other panels.

### Availability

Results in 3-4 weeks

### Gene set description

Genes in the Congenital Disorders of Glycosylation Panel and their clinical significance

Gene	Associated phenotypes	Inheritance	ClinVar	HGMD
<a href="#">ALG1*</a>	Congenital disorder of glycosylation	AR	25	43
<a href="#">ALG11*</a>	Congenital disorder of glycosylation	AR	11	14
ALG12	Congenital disorder of glycosylation	AR	11	15
ALG13	Congenital disorder of glycosylation	XL	5	12
ALG2	Congenital disorder of glycosylation, Myasthenic syndrome, congenital	AR	5	5
ALG3	Congenital disorder of glycosylation	AR	9	18

# Blueprint Genetics

ALG6	Congenital disorder of glycosylation	AR	28	24
ALG8	Congenital disorder of glycosylation	AD/AR	10	17
ALG9	Congenital disorder of glycosylation, Gillessen-Kaesbach-Nishimura syndrome	AR	4	4
ATP6V0A2	Cutis laxa, Wrinkly skin syndrome	AR	16	56
B3GLCT	Peters-plus syndrome	AR	9	15
B4GALT1	Congenital disorder of glycosylation	AR	1	2
COG1	Congenital disorder of glycosylation	AR	4	3
COG4	Congenital disorder of glycosylation	AR	12	4
COG5	Congenital disorder of glycosylation	AR	4	13
COG6	Congenital disorder of glycosylation, Shaheen syndrome	AR	10	9
COG7	Congenital disorder of glycosylation	AR	7	5
COG8	Congenital disorder of glycosylation	AR	5	7
DDOST	Congenital disorder of glycosylation	AR	3	2
DHDDS	Retinitis pigmentosa	AR	5	8
DOLK	Congenital disorder of glycosylation	AR	8	11
DPAGT1	Congenital disorder of glycosylation, Myasthenic syndrome, congenital	AR	16	32
DPM1	Congenital disorder of glycosylation	AR	9	8
DPM2	Congenital disorder of glycosylation	AR	2	2
DPM3	Congenital disorder of glycosylation, Dilated cardiomyopathy (DCM), Limb-girdle muscular dystrophy	AR	3	2
FUT8	Congenital disorder of glycosylation	AR	4	4
GMPPA	Alacrima, achalasia, and mental retardation syndrome	AR	6	12
GNE	Inclusion body myopathy, Nonaka myopathy, Sialuria	AD/AR	78	214
MAGT1	Immunodeficiency, with magnesium defect, Epstein-Barr virus infection and neoplasia, Mental retardation, X-linked 95	XL	8	14
MAN1B1	Mental retardation	AR	8	26
MGAT2	Congenital disorder of glycosylation	AR	6	5
MOGS	Congenital disorder of glycosylation	AR	7	8
MPDU1	Congenital disorder of glycosylation	AR	7	7
MPI	Congenital disorder of glycosylation	AR	27	20
NGLY1	Congenital disorder of deglycosylation	AR	26	22

PGM1	Congenital disorder of glycosylation	AR	11	35
PMM2	Congenital disorder of glycosylation	AR	76	128
RFT1	Congenital disorder of glycosylation	AR	11	13
SEC23B	Anemia, dyserythropoietic congenital	AR	18	121
SLC35A1	Congenital disorder of glycosylation	AR	4	5
SLC35A2	Congenital disorder of glycosylation	XL	16	16
SLC35C1	Congenital disorder of glycosylation, Leukocyte adhesion deficiency	AR	6	7
<u>SRD5A3*</u>	Kahrizi syndrome, Congenital disorder of glycosylation, Retinal dystrophy	AR	13	16
SSR4	Congenital disorder of glycosylation	XL	5	7
STT3A	Congenital disorder of glycosylation	AR	1	2
STT3B	Congenital disorder of glycosylation	AR	1	4
TMEM165	Congenital disorder of glycosylation	AR	4	6
TUSC3	Mental retardation	AR	6	16

\*Some regions of the gene are duplicated in the genome. [Read more.](#)

# The gene has suboptimal coverage (means <90% of the gene's target nucleotides are covered at >20x with mapping quality score (MQ>20) reads), and/or the gene has exons listed under Test limitations section that are not included in the panel as they are not sufficiently covered with high quality sequence reads.

The sensitivity to detect variants may be limited in genes marked with an asterisk (\*) or number sign (#)

Gene refers to the HGNC approved gene symbol; Inheritance refers to inheritance patterns such as autosomal dominant (AD), autosomal recessive (AR), mitochondrial (mi), X-linked (XL), X-linked dominant (XLD) and X-linked recessive (XLR); ClinVar refers to the number of variants in the gene classified as pathogenic or likely pathogenic in this database ([ClinVar](#)); HGMD refers to the number of variants with possible disease association in the gene listed in Human Gene Mutation Database ([HGMD](#)). The list of associated, gene specific phenotypes are generated from [CGD](#) or Mitomap databases.

## Non-coding disease causing variants covered by the panel

Gene	Genomic location HG19	HGVS	RefSeq	RS-number
ALG6	Chr1:63871975	c.347-13C>G	NM_013339.3	
COG5	Chr7:106898843	c.1669-15A>G	NM_006348.3	
COG6	Chr13:40273614	c.1167-24A>G	NM_020751.2	rs730882236
DHDDS	Chr1:26774026	c.441-24A>G	NM_024887.3	rs764831063
PGM1	Chr1:64113966	c.1199-222G>T	NM_001172818.1	
PMM2	Chr16:8891573		NM_000303.2	
PMM2	Chr16:8898599	c.179-25A>G	NM_000303.2	rs760689221

PMM2	Chr16:8926102	c.640-15479C>T	NM_000303.2	rs1258107584
PMM2	Chr16:8941558	c.640-23A>G	NM_000303.2	
SEC23B	Chr20:18488060	c.-571A>G	NM_006363.4	rs559854357
SEC23B	Chr20:18488615	c.-16A>G	NM_006363.4	
SEC23B	Chr20:18491731	c.221+31A>G	NM_006363.4	
SEC23B	Chr20:18491863	c.221+163A>G	NM_006363.4	rs573898514
SEC23B	Chr20:18492791	c.222-78C>T	NM_006363.4	rs150393520
SEC23B	Chr20:18526845	c.1743+168A>G	NM_006363.4	rs111951711
STT3B	Chr3:31663820	c.1539+20G>T	NM_178862.1	
TMEM165	Chr4:56284334	c.792+182G>A	NM_018475.4	rs793888506

## Test Strengths

### The strengths of this test include:

- CAP accredited laboratory
- CLIA-certified personnel performing clinical testing in a CLIA-certified laboratory
- Powerful sequencing technologies, advanced target enrichment methods and precision bioinformatics pipelines ensure superior analytical performance
- Careful construction of clinically effective and scientifically justified gene panels
- Some of the panels include the whole mitochondrial genome (please see the Panel Content section)
- Our Nucleus online portal providing transparent and easy access to quality and performance data at the patient level
- Our publicly available analytic validation demonstrating complete details of test performance
- ~2,000 non-coding disease causing variants in our clinical grade NGS assay for panels (please see 'Non-coding disease causing variants covered by this panel' in the Panel Content section)
- Our rigorous variant classification scheme
- Our systematic clinical interpretation workflow using proprietary software enabling accurate and traceable processing of NGS data
- Our comprehensive clinical statements

## Test Limitations

The following exons are not included in the panel as they are not sufficiently covered with high quality sequence reads: *ALG8* (NM\_001007027:13). Genes with suboptimal coverage in our assay are marked with number sign (#) and genes with partial, or whole gene, segmental duplications in the human genome are marked with an asterisk (\*) if they overlap with the UCSC pseudogene regions. Gene is considered to have suboptimal coverage when >90% of the gene's target nucleotides are not covered at >20x with mapping quality score (MQ>20) reads. The technology may have limited sensitivity to detect variants in genes marked with these symbols (please see the Panel content table above).

### This test does not detect the following:

- Complex inversions
- Gene conversions
- Balanced translocations
- Some of the panels include the whole mitochondrial genome but not all (please see the Panel Content section)
- Repeat expansion disorders unless specifically mentioned
- Non-coding variants deeper than  $\pm 20$  base pairs from exon-intron boundary unless otherwise indicated (please see

above Panel Content / non-coding variants covered by the panel).

## This test may not reliably detect the following:

- Low level mosaicism in nuclear genes (variant with a minor allele fraction of 14.6% is detected with 90% probability)
- Stretches of mononucleotide repeats
- Low level heteroplasmy in mtDNA (>90% are detected at 5% level)
- Indels larger than 50bp
- Single exon deletions or duplications
- Variants within pseudogene regions/duplicated segments
- Some disease causing variants present in mtDNA are not detectable from blood, thus post-mitotic tissue such as skeletal muscle may be required for establishing molecular diagnosis.

The sensitivity of this test may be reduced if DNA is extracted by a laboratory other than Blueprint Genetics.

For additional information, please refer to the Test performance section and see our Analytic Validation.

## Test performance

The Blueprint Genetics congenital disorders of glycosylation panel covers classical genes associated with disorder of glycoprotein metabolism, GEN-related myopathy and ATP6VOA2-related cutis laxa. The genes on the panel have been carefully selected based on scientific literature, mutation databases and our experience.

Our panels are sliced from our high-quality whole exome sequencing data. Please see our sequencing and detection performance table for different types of alterations at the whole exome level (Table).

Assays have been validated for different starting materials including EDTA-blood, isolated DNA (no FFPE), saliva and dry blood spots (filter card) and all provide high-quality results. The diagnostic yield varies substantially depending on the assay used, referring healthcare professional, hospital and country. Blueprint Genetics' Plus Analysis (Seq+Del/Dup) maximizes the chance to find a molecular genetic diagnosis for your patient although Sequence Analysis or Del/Dup Analysis may be a cost-effective first line test if your patient's phenotype is suggestive of a specific mutation type.

The genes on the panel have been carefully selected based on scientific literature, mutation databases and our experience.

Our panels are sectioned from our high-quality, clinical grade NGS assay. Please see our sequencing and detection performance table for details regarding our ability to detect different types of alterations (Table).

	Sensitivity % (TP/(TP+FN))	Specificity %
Single nucleotide variants	99.89% (99,153/99,266)	>99.9999%
Insertions, deletions and indels by sequence analysis		
1-10 bps	96.9% (7,563/7,806)	>99.9999%
11-50 bps	99.13% (2,524/2,546)	>99.9999%
Copy number variants (exon level dels/dups)		
1 exon level deletion (heterozygous)	100% (20/20)	NA
1 exon level deletion (homozygous)	100% (5/5)	NA
1 exon level deletion (het or homo)	100% (25/25)	NA
2-7 exon level deletion (het or homo)	100% (44/44)	NA
1-9 exon level duplication (het or homo)	75% (6/8)	NA



Simulated CNV detection		
5 exons level deletion/duplication	98.7%	100.00%
Microdeletion/-duplication sdrs (large CNVs, n=37)		
Size range (0.1-47 Mb)	100% (37/37)	

The performance presented above reached by Blueprint Genetics high-quality, clinical grade NGS sequencing assay with the following coverage metrics

Mean sequencing depth	143X
Nucleotides with >20x sequencing coverage (%)	99.86%

### Performance of Blueprint Genetics Mitochondrial Sequencing Assay.

	Sensitivity % (TP/(TP+FN))	Specificity
ANALYTIC VALIDATION (NA samples; n=4)		
Single nucleotide variants		
Heteroplasmic (45-100%)	100.0% (50/50)	100.0%
Heteroplasmic (35-45%)	100.0% (87/87)	100.0%
Heteroplasmic (25-35%)	100.0% (73/73)	100.0%
Heteroplasmic (15-25%)	100.0% (77/77)	100.0%
Heteroplasmic (10-15%)	100.0% (74/74)	100.0%
Heteroplasmic (5-10%)	100.0% (3/3)	100.0%
Heteroplasmic (<5%)	50.0% (2/4)	100.0%
CLINICAL VALIDATION (n=76 samples)		
All types		
Single nucleotide variants n=2084 SNVs		
Heteroplasmic (45-100%)	100.0% (1940/1940)	100.0%
Heteroplasmic (35-45%)	100.0% (4/4)	100.0%
Heteroplasmic (25-35%)	100.0% (3/3)	100.0%
Heteroplasmic (15-25%)	100.0% (3/3)	100.0%
Heteroplasmic (10-15%)	100.0% (9/9)	100.0%
Heteroplasmic (5-10%)	92.3%(12/13)	99.98%
Heteroplasmic (<5%)	88.7% (47/53)	99.93%



# Blueprint Genetics



Insertions and deletions by sequence analysis n=42 indels		
Heteroplasmic (45-100%) 1-10bp	100.0% (32/32)	100.0%
Heteroplasmic (5-45%) 1-10bp	100.0% (3/3)	100.0%
Heteroplasmic (<5%) 1-10bp	100.0% (5/5)	>0.9999
SIMULATION DATA /(mitomap mutations)		
Insertions, and deletions 1-24 bps by sequence analysis; n=17		
Homoplasmic (100%) 1-24bp	100.0% (17/17)	99.98%
Heteroplasmic (50%)	100.0% (17/17)	99.99%
Heteroplasmic (25%)	100.0% (17/17)	100.0%
Heteroplasmic (20%)	100.0% (17/17)	100.0%
Heteroplasmic (15%)	100.0% (17/17)	100.0%
Heteroplasmic (10%)	94.1% (16/17)	100.0%
Heteroplasmic (5%)	94.1% (16/17)	100.0%
Copy number variants (separate artificial mutations; n=1500)		
Homoplasmic (100%) 500 bp, 1kb, 5 kb	100.0%	100.0%
Heteroplasmic (50%) 500 bp, 1kb, 5 kb	100.0%	100.0%
Heteroplasmic (30%) 500 bp, 1kb, 5 kb	100.0%	100.0%
Heteroplasmic (20%) 500 bp, 1kb, 5 kb	99.7%	100.0%
Heteroplasmic (10%) 500 bp, 1kb, 5 kb	99.0%	100.0%
The performance presented above reached by following coverage metrics at assay level (n=66)		
	Mean of medians	Median of medians
Mean sequencing depth MQ0 (clinical)	18224X	17366X
Nucleotides with >1000x MQ0 sequencing coverage (%) (clinical)	100%	
rho zero cell line (=no mtDNA), mean sequencing depth	12X	

## Bioinformatics

The target region for each gene includes coding exons and ±20 base pairs from the exon-intron boundary. In addition, the panel includes non-coding variants if listed above (Non-coding variants covered by the panel). Some regions of the gene(s) may be removed from the panel if specifically mentioned in the ‘Test limitations’ section above. The sequencing data generated in our laboratory is analyzed with our proprietary data analysis and annotation pipeline, integrating state-of-the-art algorithms and industry-standard software solutions. Incorporation of rigorous quality control steps throughout the workflow of the pipeline ensures the consistency, validity and accuracy of results. Our pipeline is streamlined to maximize sensitivity without sacrificing specificity. We have incorporated a number of reference population databases and mutation databases such as, but not limited, to [1000 Genomes Project](#), [gnomAD](#), [ClinVar](#) and [HGMD](#) into our clinical interpretation software to



# Blueprint Genetics

make the process effective and efficient. For missense variants, *in silico* variant prediction tools such as SIFT, PolyPhen, MutationTaster are used to assist with variant classification. Through our online ordering and statement reporting system, Nucleus, the customer has an access to details of the analysis, including patient specific sequencing metrics, a gene level coverage plot and a list of regions with inadequate coverage if present. This reflects our mission to build fully transparent diagnostics where customers have easy access to crucial details of the analysis process.

## Clinical interpretation

We provide customers with the most comprehensive clinical report available on the market. Clinical interpretation requires a fundamental understanding of clinical genetics and genetic principles. At Blueprint Genetics, our PhD molecular geneticists, medical geneticists and clinical consultants prepare the clinical statement together by evaluating the identified variants in the context of the phenotypic information provided in the requisition form. Our goal is to provide clinically meaningful statements that are understandable for all medical professionals regardless of whether they have formal training in genetics.

Variant classification is the corner stone of clinical interpretation and resulting patient management decisions. Our classifications follow the [Blueprint Genetics Variant Classification Schemes](#) based on the [ACMG guideline 2015](#). Minor modifications were made to increase reproducibility of the variant classification and improve the clinical validity of the report. Our experience with tens of thousands of clinical cases analyzed at our laboratory allowed us to further develop the industry standard.

The final step in the analysis of sequence variants is confirmation of variants classified as pathogenic or likely pathogenic using bi-directional Sanger sequencing. Variant(s) fulfilling the following criteria are not Sanger confirmed: the variant quality score is above the internal threshold for a true positive call, and visual check-up of the variant at IGV is in-line with the variant call. Reported variants of uncertain significance are confirmed with bi-directional Sanger sequencing only if the quality score is below our internally defined quality score for true positive call. Reported copy number variations with a size <10 exons are confirmed by orthogonal methods such as qPCR if the specific CNV has been seen less than three times at Blueprint Genetics.

Our clinical statement includes tables for sequencing and copy number variants that include basic variant information (genomic coordinates, HGVS nomenclature, zygosity, allele frequencies, in silico predictions, OMIM phenotypes and classification of the variant). In addition, the statement includes detailed descriptions of the variant, gene and phenotype(s) including the role of the specific gene in human disease, the mutation profile, information about the gene's variation in population cohorts and detailed information about related phenotypes. We also provide links to the references used, congress abstracts and mutation variant databases used to help our customers further evaluate the reported findings if desired. The conclusion summarizes all of the existing information and provides our rationale for the classification of the variant.

Identification of pathogenic or likely pathogenic variants in dominant disorders or their combinations in different alleles in recessive disorders are considered molecular confirmation of the clinical diagnosis. In these cases, family member testing can be used for risk stratification within the family. In the case of variants of uncertain significance (VUS), we do not recommend family member risk stratification based on the VUS result. Furthermore, in the case of VUS, we do not recommend the use of genetic information in patient management or genetic counseling.

Our interpretation team analyzes millions of variants from thousands of individuals with rare diseases. Thus, our database, and our understanding of variants and related phenotypes, is growing by leaps and bounds. Our laboratory is therefore well positioned to re-classify previously reported variants as new information becomes available. If a variant previously reported by Blueprint Genetics is re-classified, our laboratory will issue a follow-up statement to the original ordering health care provider at no additional cost.

#}

## ICD codes

Commonly used ICD-10 codes when ordering the Congenital Disorders of Glycosylation Panel

ICD-10	Disease
E77.9	Disorder of glycoprotein metabolism



## Accepted sample types

- EDTA blood, min. 1 ml
- Purified DNA, min. 3µg\*
- Saliva (Oragene DNA OG-500 kit)

Label the sample tube with your patient's name, date of birth and the date of sample collection.

Note that we do not accept DNA samples isolated from formalin-fixed paraffin-embedded (FFPE) tissue.

## Resources

- [CDG Care](#)
- [Foundation Glycosylation](#)
- [GeneReviews - ATP6V0A2-Related Cutis Laxa](#)
- [GeneReviews - Congenital Disorders of Glycosylation](#)
- [GeneReviews - Congenital Disorders of N-Linked Glycosylation and Multiple Pathway Overview.](#)
- [GeneReviews- \\*ATP6V0A2\\*-Related Cutis Laxa](#)
- [National Organization for Rare Disorders](#)
- [Portuguese Association CDG](#)