

Dementia Panel

Test code: NE2301

Is a 58 gene panel that includes assessment of non-coding variants.

In addition, it also includes the maternally inherited mitochondrial genome.

Is ideal for patients with a clinical suspicion of dementia.

About Dementia

Frontotemporal dementia (FTD) is a clinically and pathologically heterogeneous group of non-Alzheimer dementias characterized by selective, progressive cortical atrophy involving the frontal or temporal lobes. FTD is substantially less common than Alzheimer's disease, with estimates of population prevalence ranging from 4-5 per 100,000 before age 65. Age of onset is typically in the sixth decade of life. However, it may begin as early 30 or as late as the ninth decade. Approximately 20-50% of individuals with FTD have an affected first degree relative. FTD has a substantial genetic component, with an autosomal dominant or X-linked inheritance pattern. It is estimated that 10% of patients with FTD have a disease causing mutation in a single gene. The *APOE* E4 haplotype confers a significant risk for Alzheimer's disease related dementia, especially in homozygous state. Therefore, this haplotype is reported from this panel, if detected in homozygous state.

Availability

4 weeks

Gene Set Description

Genes in the Dementia Panel and their clinical significance

Gene	Associated phenotypes	Inheritance	ClinVar	HGMD
ABCA7	Alzheimer disease	AD	1	139
APOE	Sea-blue histiocyte disease, Dysbetalipoproteinemia, familial (Hyperlipoproteinemia), Lipoprotein glomerulopathy	AD/AR	31	55
APP	Alzheimer disease, Cerebral amyloid angiopathy	AD/AR	18	100
CHMP2B	Amyotrophic lateral sclerosis, CHMP2B-related, Frontotemporal dementia	AD	6	21
CSF1R	Leukoencephalopathy, diffuse hereditary, with spheroids	AD	56	83
FUS	Amyotrophic lateral sclerosis, Essential tremor	AD/AR	22	111
GRN	Frontotemporal lobar degeneration with TDP43 inclusions, GRN-related, Neuronal ceroid lipofuscinosis	AD/AR	43	214
MAPT	Pick disease, Frontotemporal dementia, Parkinson-dementia syndrome, Supranuclear palsy, progressive	AD/AR	26	104
MT-ATP6	Neuropathy, ataxia, and retinitis pigmentosa, Leber hereditary optic neuropathy, Ataxia and polyneuropathy, adult-onset, Cardiomyopathy, infantile hypertrophic, Leigh syndrome, Striatonigral degeneration, infantile, mitochondrial	Mitochondrial	19	
MT-ATP8	Cardiomyopathy, apical hypertrophic, and neuropathy, Cardiomyopathy, infantile hypertrophic	Mitochondrial	4	

Blueprint Genetics

MT-CO1	Myoglobinuria, recurrent, Leber hereditary optic neuropathy, Sideroblastic anemia, Cytochrome C oxidase deficiency, Deafness, mitochondrial	Mitochondrial	17
MT-CO2	Cytochrome c oxidase deficiency	Mitochondrial	8
MT-CO3	Cytochrome c oxidase deficiency, Leber hereditary optic neuropathy	Mitochondrial	9
MT-CYB		Mitochondrial	69
MT-ND1	Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes, Leber hereditary optic neuropathy, Leber optic atrophy and dystonia	Mitochondrial	21
MT-ND2	Leber hereditary optic neuropathy, Mitochondrial complex I deficiency	Mitochondrial	6
MT-ND3	Leber optic atrophy and dystonia, Mitochondrial complex I deficiency	Mitochondrial	7
MT-ND4	Leber hereditary optic neuropathy, Leber optic atrophy and dystonia, Mitochondrial complex I deficiency	Mitochondrial	11
MT-ND4L	Leber hereditary optic neuropathy	Mitochondrial	2
MT-ND5	Myoclonic epilepsy with ragged red fibers, Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes, Leber hereditary optic neuropathy, Mitochondrial complex I deficiency	Mitochondrial	19
MT-ND6	Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes, Oncocytoma, Leber hereditary optic neuropathy, Leber optic atrophy and dystonia, Mitochondrial complex I deficiency	Mitochondrial	16
MT-RNR1	Deafness, mitochondrial	Mitochondrial	3
MT-RNR2	Chloramphenicol toxicity/resistance	Mitochondrial	2
MT-TA		Mitochondrial	4
MT-TC	Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes	Mitochondrial	3
MT-TD		Mitochondrial	1
MT-TE	Diabetes-deafness syndrome, Mitochondrial myopathy, infantile, transient, Mitochondrial myopathy with diabetes	Mitochondrial	5
MT-TF	Myoclonic epilepsy with ragged red fibers, Nephropathy, tubulointerstitial, Encephalopathy, mitochondrial, Epilepsy, mitochondrial, Myopathy, mitochondrial, Mitochondrial encephalomyopathy with lactic acidosis and stroke-like episodes	Mitochondrial	7
MT-TG		Mitochondrial	3
MT-TH		Mitochondrial	4
MT-TI		Mitochondrial	7
MT-TK	Myoclonic epilepsy with ragged red fibers, Leigh syndrome	Mitochondrial	5
MT-TL1	Cytochrome c oxidase deficiency, Myoclonic epilepsy with ragged red fibers, Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes, Diabetes-deafness syndrome, Cyclic vomiting syndrome, SIDS, susceptibility to	Mitochondrial	14

Blueprint Genetics

MT-TL2	Mitochondrial multisystemic disorder, Progressive external ophthalmoplegia, Mitochondrial Myopathy, Mitochondrial encephalomyopathy with lactic acidosis and stroke-like episodes	Mitochondrial	5	
MT-TM	Leigh syndrome, Mitochondrial multisystemic disorder	Mitochondrial	1	
MT-TN	Progressive external ophthalmoplegia, Mitochondrial multisystemic disorder	Mitochondrial	3	
MT-TP		Mitochondrial	2	
MT-TQ	Mitochondrial multisystemic disorder	Mitochondrial	2	
MT-TR	Encephalopathy, mitochondrial	Mitochondrial	2	
MT-TS1	Myoclonic epilepsy with ragged red fibers, Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes	Mitochondrial	10	
MT-TS2	Mitochondrial multisystemic disorder	Mitochondrial	2	
MT-TT		Mitochondrial	5	
MT-TV	Hypertrophic cardiomyopathy (HCM), Leigh syndrome, Mitochondrial multisystemic disorder, Mitochondrial encephalomyopathy with lactic acidosis and stroke-like episodes	Mitochondrial	3	
MT-TW	Leigh syndrome, Myopathy, mitochondrial	Mitochondrial	8	
MT-TY	Mitochondrial multisystemic disorder	Mitochondrial	4	
PRNP	Dementia, Lewy body, Creutzfeldt-Jakob disease, Huntington disease-like, Gerstmann-Straussler disease, Spongiform encephalopathy with neuropsychiatric features, Insomnia, fatal familial	AD/AR	26	104
PSEN1	Dilated cardiomyopathy (DCM), Acne inversa, familial, 3, Dementia, frontotemporal, Pick disease, Alzheimer disease	AD	57	306
PSEN2	Peripartum/pregnancy-associated cardiomyopathy, Dilated cardiomyopathy (DCM), Alzheimer disease, 4	AD	9	60
RNF216*	Cerebellar ataxia and hypogonadotropic hypogonadism (Gordon Holmes syndrome)	AR	10	14
SIGMAR1	Amyotrophic lateral sclerosis, Spinal muscular atrophy, distal, Frontotemporal lobar degeneration-motor neuron disease	AR	6	14
SNCA	Parkinson disease, Dementia with Lewy bodies	AD	7	35
SORL1	Early-onset Alzheimer disease	AD	3	134
TARDBP*	Amyotrophic lateral sclerosis	AD	20	69
TREM2	Nasu-Hakola disease, Early-onset dementia without bone cysts, Polycystic lipomembranous osteodysplasia with sclerosing leukoencephalopathy	AR	14	48
TUBA4A	Amyotrophic lateral sclerosis 22	AD	6	13
UBE3A*	Angelman syndrome	AD	176	202
UBQLN2	Amyotrophic lateral sclerosis	XL	5	31

VCP	Amyotrophic lateral sclerosis, Inclusion body myopathy with early-onset Paget disease, Charcot-Marie-Tooth disease	AD	17	61
-----	--	----	----	----

*Some regions of the gene are duplicated in the genome. [Read more.](#)

The gene has suboptimal coverage (means <90% of the gene's target nucleotides are covered at >20x with mapping quality score (MQ>20) reads), and/or the gene has exons listed under Test limitations section that are not included in the panel as they are not sufficiently covered with high quality sequence reads.

The sensitivity to detect variants may be limited in genes marked with an asterisk (*) or number sign (#). Due to possible limitations these genes may not be available as single gene tests.

Gene refers to the HGNC approved gene symbol; Inheritance refers to inheritance patterns such as autosomal dominant (AD), autosomal recessive (AR), mitochondrial (mi), X-linked (XL), X-linked dominant (XLD) and X-linked recessive (XLR); ClinVar refers to the number of variants in the gene classified as pathogenic or likely pathogenic in this database ([ClinVar](#)); HGMD refers to the number of variants with possible disease association in the gene listed in Human Gene Mutation Database ([HGMD](#)). The list of associated, gene specific phenotypes are generated from [CGD](#) or Mitomap databases.

Non-coding disease causing variants covered by the panel

Gene	Genomic location HG19	HGVS	RefSeq	RS-number
APP	Chr21:27253648	c.*331_*332delAT	NM_000484.3	
APP	Chr21:27543203	c.-49+100C>A	NM_001136131.2	
APP	Chr21:27543454	c.-516C>G	NM_000484.3	rs539645405
APP	Chr21:27543619	c.-681G>A	NM_000484.3	rs187510057
CSF1R	Chr5:149440654	c.1859-119G>A	NM_005211.3	
FUS	Chr16:31202807	c.*48G>A	NM_004960.3	rs376510148
FUS	Chr16:31202818	c.*59G>A	NM_004960.3	
FUS	Chr16:31202863	c.*105dupT	NM_004960.3	
FUS	Chr16:31202867	c.*108C>T	NM_004960.3	rs780606789
FUS	Chr16:31202869	c.*110G>A	NM_004960.3	
FUS	Chr16:31202891	c.*132C>A	NM_004960.3	rs565540429
FUS	Chr16:31202949	c.*190C>A	NM_004960.3	
GRN	Chr17:42422701	c.-9A>G	NM_002087.2	
GRN	Chr17:42422705	c.-8+3A>T	NM_002087.2	rs63751020
GRN	Chr17:42422705	c.-8+3A>G	NM_002087.2	
GRN	Chr17:42422707	c.-8+5G>C	NM_002087.2	rs63750313
MAPT	Chr17:44087661	c.1774-15T>C	NM_016835.4	
MAPT	Chr17:44087779	c.1866+11T>C	NM_016835.4	rs63751394

Blueprint Genetics

MAPT	Chr17:44087780	c.1866+12C>T	NM_016835.4	rs63750916
MAPT	Chr17:44087781	c.1866+13A>G	NM_016835.4	rs63750308
MAPT	Chr17:44087782	c.1866+14C>T	NM_016835.4	rs63750972
MAPT	Chr17:44087783	c.1866+15A>C	NM_016835.4	
MAPT	Chr17:44087784	c.1866+16C>T	NM_016835.4	rs63751011
MAPT	Chr17:44087787	c.1866+19C>G	NM_016835.4	rs63750162
PSEN1	Chr14:73673071	c.869-23_869-22insTGGAATTTTGTGCTGTTG	NM_000021.3	
SIGMAR1	Chr9:34635578	c.*51G>T	NM_005866.2	rs768783740
SNCA	Chr4:90647315	c.*464C>A	NM_000345.3	rs183204610
TARDBP	Chr1:11082794	c.*83T>C	NM_007375.3	rs80356744
TARDBP	Chr1:11083408	c.*697G>A	NM_007375.3	rs387906334
VCP	Chr9:35072710	c.-360G>C	NM_007126.3	

Test Strengths

The *APOE* E4/E4 haplotype is reported if identified on this panel.

The strengths of this test include:

- CAP accredited laboratory
- CLIA-certified personnel performing clinical testing in a CLIA-certified laboratory
- Powerful sequencing technologies, advanced target enrichment methods and precision bioinformatics pipelines ensure superior analytical performance
- Careful construction of clinically effective and scientifically justified gene panels
- Some of the panels include the whole mitochondrial genome (please see the Panel Content section)
- Our Nucleus online portal providing transparent and easy access to quality and performance data at the patient level
- Our publicly available analytic validation demonstrating complete details of test performance
- ~2,000 non-coding disease causing variants in our clinical grade NGS assay for panels (please see 'Non-coding disease causing variants covered by this panel' in the Panel Content section)
- Our rigorous variant classification scheme
- Our systematic clinical interpretation workflow using proprietary software enabling accurate and traceable processing of NGS data
- Our comprehensive clinical statements

Test Limitations

Genes with partial, or whole gene, segmental duplications in the human genome are marked with an asterisk (*) if they overlap with the UCSC pseudogene regions. The technology may have limited sensitivity to detect variants in genes marked with these symbols (please see the Panel content table above).

This test does not detect the following:

- Complex inversions

Blueprint Genetics

- Gene conversions
- Balanced translocations
- Some of the panels include the whole mitochondrial genome (please see the Panel Content section)
- Repeat expansion disorders unless specifically mentioned
- Non-coding variants deeper than ± 20 base pairs from exon-intron boundary unless otherwise indicated (please see above Panel Content / non-coding variants covered by the panel).

This test may not reliably detect the following:

- Low level mosaicism in nuclear genes (variant with a minor allele fraction of 14.6% is detected with 90% probability)
- Stretches of mononucleotide repeats
- Low level heteroplasmy in mtDNA (>90% are detected at 5% level)
- Indels larger than 50bp
- Single exon deletions or duplications
- Variants within pseudogene regions/duplicated segments
- Some disease causing variants present in mtDNA are not detectable from blood, thus post-mitotic tissue such as skeletal muscle may be required for establishing molecular diagnosis.

The sensitivity of this test may be reduced if DNA is extracted by a laboratory other than Blueprint Genetics.

For additional information, please refer to the Test performance section and see our Analytic Validation.

Test Performance

The genes on the panel have been carefully selected based on scientific literature, mutation databases and our experience.

Our panels are sectioned from our high-quality, clinical grade NGS assay. Please see our sequencing and detection performance table for details regarding our ability to detect different types of alterations (Table).

Assays have been validated for various sample types including EDTA-blood, isolated DNA (excluding from formalin fixed paraffin embedded tissue), saliva and dry blood spots (filter cards). These sample types were selected in order to maximize the likelihood for high-quality DNA yield. The diagnostic yield varies depending on the assay used, referring healthcare professional, hospital and country. Plus analysis increases the likelihood of finding a genetic diagnosis for your patient, as large deletions and duplications cannot be detected using sequence analysis alone. Blueprint Genetics' Plus Analysis is a combination of both sequencing and deletion/duplication (copy number variant (CNV)) analysis.

The performance metrics listed below are from an initial validation performed at our main laboratory in Finland. The performance metrics of our laboratory in Seattle, WA, are equivalent.

Performance of Blueprint Genetics high-quality, clinical grade NGS sequencing assay for panels.

	Sensitivity % (TP/(TP+FN))	Specificity %
Single nucleotide variants	99.89% (99,153/99,266)	>99.9999%
Insertions, deletions and indels by sequence analysis		
1-10 bps	99.2% (7,745/7,806)	>99.9999%
11-50 bps	99.13% (2,524/2,546)	>99.9999%
Copy number variants (exon level dels/dups)		
1 exon level deletion (heterozygous)	100% (20/20)	NA
1 exon level deletion (homozygous)	100% (5/5)	NA
1 exon level deletion (het or homo)	100% (25/25)	NA

Blueprint Genetics



2-7 exon level deletion (het or homo)	100% (44/44)	NA
1-9 exon level duplication (het or homo)	75% (6/8)	NA
Simulated CNV detection		
5 exons level deletion/duplication	98.7%	100.00%
Microdeletion/-duplication sdrs (large CNVs, n=37)		
Size range (0.1-47 Mb)	100% (25/25)	

The performance presented above reached by Blueprint Genetics high-quality, clinical grade NGS sequencing assay with the following coverage metrics

Mean sequencing depth	143X
Nucleotides with >20x sequencing coverage (%)	99.86%

Performance of Blueprint Genetics Mitochondrial Sequencing Assay.

	Sensitivity %	Specificity %
ANALYTIC VALIDATION (NA samples; n=4)		
Single nucleotide variants		
Heteroplasmic (45-100%)	100.0% (50/50)	100.0%
Heteroplasmic (35-45%)	100.0% (87/87)	100.0%
Heteroplasmic (25-35%)	100.0% (73/73)	100.0%
Heteroplasmic (15-25%)	100.0% (77/77)	100.0%
Heteroplasmic (10-15%)	100.0% (74/74)	100.0%
Heteroplasmic (5-10%)	100.0% (3/3)	100.0%
Heteroplasmic (<5%)	50.0% (2/4)	100.0%
CLINICAL VALIDATION (n=76 samples)		
All types		
Single nucleotide variants n=2026 SNVs		
Heteroplasmic (45-100%)	100.0% (1940/1940)	100.0%
Heteroplasmic (35-45%)	100.0% (4/4)	100.0%
Heteroplasmic (25-35%)	100.0% (3/3)	100.0%
Heteroplasmic (15-25%)	100.0% (3/3)	100.0%
Heteroplasmic (10-15%)	100.0% (9/9)	100.0%



Blueprint Genetics



Heteroplasmic (5-10%)	92.3% (12/13)	99.98%
Heteroplasmic (<5%)	88.9% (48/54)	99.93%
Insertions and deletions by sequence analysis n=40 indels		
Heteroplasmic (45-100%) 1-10bp	100.0% (32/32)	100.0%
Heteroplasmic (5-45%) 1-10bp	100.0% (3/3)	100.0%
Heteroplasmic (<5%) 1-10bp	100.0% (5/5)	99,997%
SIMULATION DATA /(mitomap mutations)		
Insertions, and deletions 1-24 bps by sequence analysis; n=17		
Homoplasmic (100%) 1-24bp	100.0% (17/17)	99.98%
Heteroplasmic (50%)	100.0% (17/17)	99.99%
Heteroplasmic (25%)	100.0% (17/17)	100.0%
Heteroplasmic (20%)	100.0% (17/17)	100.0%
Heteroplasmic (15%)	100.0% (17/17)	100.0%
Heteroplasmic (10%)	94.1% (16/17)	100.0%
Heteroplasmic (5%)	94.1% (16/17)	100.0%
Copy number variants (separate artificial mutations; n=1500)		
Homoplasmic (100%) 500 bp, 1kb, 5 kb	100.0%	100.0%
Heteroplasmic (50%) 500 bp, 1kb, 5 kb	100.0%	100.0%
Heteroplasmic (30%) 500 bp, 1kb, 5 kb	100.0%	100.0%
Heteroplasmic (20%) 500 bp, 1kb, 5 kb	99.7%	100.0%
Heteroplasmic (10%) 500 bp, 1kb, 5 kb	99.0%	100.0%
The performance presented above reached by following coverage metrics at assay level (n=66)		
	Mean of medians	Median of medians
Mean sequencing depth MQ0 (clinical)	18224X	17366X
Nucleotides with >1000x MQ0 sequencing coverage (%) (clinical)	100%	
rho zero cell line (=no mtDNA), mean sequencing depth	12X	

Bioinformatics

The target region for each gene includes coding exons and ± 20 base pairs from the exon-intron boundary. In addition, the panel includes non-coding variants if listed above (Non-coding variants covered by the panel). Some regions of the gene(s) may be removed from the panel if specifically mentioned in the 'Test limitations' section above. The sequencing data generated in our laboratory is analyzed with our proprietary data analysis and annotation pipeline, integrating state-of-the-art



Blueprint Genetics



algorithms and industry-standard software solutions. Incorporation of rigorous quality control steps throughout the workflow of the pipeline ensures the consistency, validity and accuracy of results. Our pipeline is streamlined to maximize sensitivity without sacrificing specificity. We have incorporated a number of reference population databases and mutation databases such as, but not limited to, [1000 Genomes Project](#), [gnomAD](#), [ClinVar](#) and [HGMD](#) into our clinical interpretation software to make the process effective and efficient. For missense variants, *in silico* variant prediction tools such as [SIFT](#), [PolyPhen](#), [MutationTaster](#) are used to assist with variant classification. Through our online ordering and statement reporting system, Nucleus, the customer has an access to details of the analysis, including patient specific sequencing metrics, a gene level coverage plot and a list of regions with inadequate coverage if present. This reflects our mission to build fully transparent diagnostics where customers have easy access to crucial details of the analysis process.

Clinical Interpretation

We provide customers with the most comprehensive clinical report available on the market. Clinical interpretation requires a fundamental understanding of clinical genetics and genetic principles. At Blueprint Genetics, our PhD molecular geneticists, medical geneticists and clinical consultants prepare the clinical statement together by evaluating the identified variants in the context of the phenotypic information provided in the requisition form. Our goal is to provide clinically meaningful statements that are understandable for all medical professionals regardless of whether they have formal training in genetics.

Variant classification is the corner stone of clinical interpretation and resulting patient management decisions. Our classifications follow the [ACMG guideline 2015](#).

The final step in the analysis of sequence variants is confirmation of variants classified as pathogenic or likely pathogenic using bi-directional Sanger sequencing. Variant(s) fulfilling the following criteria are not Sanger confirmed: the variant quality score is above the internal threshold for a true positive call, and visual check-up of the variant at IGV is in-line with the variant call. Reported variants of uncertain significance are confirmed with bi-directional Sanger sequencing only if the quality score is below our internally defined quality score for true positive call. Reported copy number variations with a size <10 exons are confirmed by orthogonal methods such as qPCR if the specific CNV has been seen less than three times at Blueprint Genetics.

Our clinical statement includes tables for sequencing and copy number variants that include basic variant information (genomic coordinates, HGVS nomenclature, zygosity, allele frequencies, *in silico* predictions, OMIM phenotypes and classification of the variant). In addition, the statement includes detailed descriptions of the variant, gene and phenotype(s) including the role of the specific gene in human disease, the mutation profile, information about the gene's variation in population cohorts and detailed information about related phenotypes. We also provide links to the references used, congress abstracts and mutation variant databases used to help our customers further evaluate the reported findings if desired. The conclusion summarizes all of the existing information and provides our rationale for the classification of the variant.

Identification of pathogenic or likely pathogenic variants in dominant disorders or their combinations in different alleles in recessive disorders are considered molecular confirmation of the clinical diagnosis. In these cases, family member testing can be used for risk stratification within the family. In the case of variants of uncertain significance (VUS), we do not recommend family member risk stratification based on the VUS result. Furthermore, in the case of VUS, we do not recommend the use of genetic information in patient management or genetic counseling.

Our interpretation team analyzes millions of variants from thousands of individuals with rare diseases. Thus, our database, and our understanding of variants and related phenotypes, is growing by leaps and bounds. Our laboratory is therefore well positioned to re-classify previously reported variants as new information becomes available. If a variant previously reported by Blueprint Genetics is re-classified, our laboratory will issue a follow-up statement to the original ordering health care provider at no additional cost.

CPT code(s) *

81406 x6, 81405 x2, 81404, 81479, 81460, 81465

* The CPT codes provided are based on AMA guidelines and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.

ICD Codes

Refer to the most current version of ICD-10-CM manual for a complete list of ICD-10 codes.

Sample Requirements

- Blood (min. 1ml) in an EDTA tube
- Extracted DNA, min. 2 µg in TE buffer or equivalent
- Saliva (Please see [Sample Requirements](#) for accepted saliva kits)

Label the sample tube with your patient's name, date of birth and the date of sample collection.

We do not accept DNA samples isolated from formalin-fixed paraffin-embedded (FFPE) tissue. In addition, if the patient is affected with a hematological malignancy, DNA extracted from a non-hematological source (e.g. skin fibroblasts) is strongly recommended.

Please note that, in rare cases, mitochondrial genome (mtDNA) variants may not be detectable in blood or saliva in which case DNA extracted from post-mitotic tissue such as skeletal muscle may be a better option.

Read more about our sample requirements [here](#).

For Patients

Other

- [Association for Frontotemporal Degeneration](#)
- [Dementia](#)
- [Dementia Advocacy and Support Network](#)
- [Dyer SM et al. Clinical practice guidelines and principles of care for people with dementia in Australia. Aust Fam Physician. 2016 Dec;45\(12\):884-889.](#)
- [Frontotemporal Dementia Support Group](#)
- [Frontotemporal Lobar Degeneration Association](#)
- [Rosness TA et al. Frontotemporal Dementia: An Updated Clinician's Guide. J Geriatr Psychiatry Neurol. 2016 Sep;29\(5\):271-80.](#)